



# QUICK START GUIDE

# Table of Contents


Web Portal System Setup	1
Rostering Athletes in Web Portal	2
Team App Setup	3
Assigning IMM's	4
Charging Monitors	5
Daily Usage	6
Post Session Analysis	7

# WEB PORTAL SETUP

To get started email [support@preventbiometrics.com](mailto:support@preventbiometrics.com) the email address that you wish to create an account with. Prevent staff will confirm that your email has been assigned in the system


## **LOGGING IN**

Once Prevent staff have notified you that your email has been assigned in the portal, proceed to [www.portal.preventbiometrics.com](http://www.portal.preventbiometrics.com) to create your Prevent account



The screenshot shows the registration page of the Prevent web portal. At the top, the Prevent logo and tagline "REAL TIME HEAD IMPACT MONITORING" are displayed. Below this, a white box contains the heading "PLEASE PROVIDE THE FOLLOWING DETAILS." and several input fields: "EMAIL ADDRESS" (with a "SEND VERIFICATION CODE" button), "NEW PASSWORD", "CONFIRM NEW PASSWORD", "FIRST NAME", "LAST NAME", "ORGANIZATION", and "COUNTRY/REGION" (with a dropdown arrow).

Fill Out Account Details



The screenshot shows the login page of the Prevent web portal. At the top, the Prevent logo and tagline "REAL TIME HEAD IMPACT MONITORING" are displayed. Below this, a white box contains the heading "SIGN IN WITH YOUR EXISTING ACCOUNT" and two input fields: "EMAIL ADDRESS" and "PASSWORD" (with a "FORGOT YOUR PASSWORD?" link). At the bottom of the box, there are two buttons: "SIGN IN" and "SIGN UP NOW". A yellow box highlights the "SIGN UP NOW" button, and a green arrow points from the text "Fill Out Account Details" to it.


# ROSTERING ATHLETES

The screenshot shows the Prevent portal interface for a 'Demo Team'. The left sidebar contains navigation options: Administration, All Teams, CURRENT TEAMS (expanded to show Demo Team), and Demo Team. The main content area displays team statistics: 0 IMPACTS, 0 UNSYNCD, and Hours. Below this is a table with columns for SPORT (Football), ATHLETES (0), and SEASON (Mar 1, 2021 - Jan 1, 2022). A navigation bar includes HEAD IMPACTS, COLLISIONS, ATHLETES (selected), SYSTEM STATUS, SCHEDULE, and EDIT TEAM... The 'TEAM ATHLETES' section features tabs for IMPACTS and COLLISIONS, filters for 'All season' and 'Positions: All positions', and a blue 'EDIT ROSTER' button. Below the button is a table header with columns: #, PLAYER NAME, IMPACTS, IMPACT DAYS, IMPACTS PER DAY, NO DATE IMPACTS, MED+, HIGH, and IMPACT COUNTS.


Select the “Edit Roster” tab and select the “Add Athletes” icon to begin adding athletes to your roster. The athletes added in the portal will automatically populate in the iOS Prevent team app

# TEAM APP SETUP

Using an iOS device, go to the App Store and search "Prevent Biometrics". Proceed to download the Prevent app



Upon opening the Team App, you will be presented with the following home page. "Sign In" to your account **using the same email address and password** used in the portal setup process





# ASSIGNING IMMS

Select the team you would like to view from the drop-down menu at the top of the app.

Each monitor has a unique device serial number. Prior to use, each monitor must be assigned to a player

Select the "+ Assign Monitor" icon next to the athlete you wish to assign a monitor to. Scan the serial number barcode to finalize the monitor assignment



# CHARGING MONITORS

Place the IMM in the charging case and securely close the lid. Ensure the IMM is pushed all the way down in the charging case, and the back left section of the IMM is touching wall of the solo case. Ensure each solo case is labeled with a key identifier to each athlete



Plug the charging case into any slot in the team case. Refrain from jamming the case into the slot with force. Ensure the team case is plugged into a power source



Use the iOS app to confirm that all monitors have fully charged and are connecting to the iOS app prior to being deployed into the field



## INDICATOR LIGHTS

- CHARGING
- IN CASE BUT NOT ALIGNED, INVESTIGATE
- NOT IN CASE, OR MALFUNCTIONING

Use iOS app to confirm monitor functionality for Blue & Red Solo Cases

Please contact [support@preventbiometrics.com](mailto:support@preventbiometrics.com) with any questions

# DAILY USAGE

Ensure all monitors have charged for at least 70 minutes. Ensure proper alignment in solo case using the solo case light indicators



Open the Prevent iOS app while within range of the charging monitors for 1 minute prior to use



Athletes remove their monitor from solo case and wear during session




(Optional) Open iOS app during session: IMM's within range (300 feet) of the iOS app will sync impact data in real-time



Following each session, ensure all devices are returned to their charging position



Following the session, open the Prevent iOS app while within range of the charging monitors for 10 minutes to sync any remaining data





# POST SESSION ANALYSIS

Analyze session data through a wide-range of reporting features in the Prevent portal and iOS application.

Automated reports are sent to your email inbox at 6 AM every morning following a data collection session

